EXERGEN

C O R P O R A T I O N

CLINICAL STUDIES OVERVIEW

Exergen Corporation
400 Pleasant Street
Watertown, MA 02472
USA

Email: service@exergen.com
Phone: 617.923.9900

[< Back to categories](#)

CL-901-CH-V0

CATEGORIES

Adult

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Adult CCU

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School Nurse

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White Papers

White Paper #1: Measuring body temperature:
A simple, effective and reliable solution

Case Studies

Case Study #1: Hospitals Save 1,4 million Euro By Adopting TAT- 5000 Temporal Artery Thermometer
Case Study #2: Appeals Court Upholds Exergen Temporal Artery Thermometer Patent Infringement by Helen of Troy's Braun and Vicks Forehead Thermometers
Case Study #3: Temporal Artery Thermometry very comfortable for preterm neonates
Case Study #4: University Medical Centre Groningen chooses TAT-5000 of Exergen Medical as standard thermometer throughout hospital
Case Study #5: 90 percent of nurses prefer thermometer connected to vital signs monitor
Case Study #6: Why Seattle Children's Hospital decided to adopt temporal artery thermometers
Case Study #7: UL grants Exergen Corporation ISO 13485:2016 certificate
Case Study #8: After UMCG, Beatrixoord, Centre for Rehabilitation in Haren, also chooses infrared thermometers from Exergen
Case Study #9: Exergen medical forms strategic distribution partnership with international medical trade consultants of Israel
Case Study #10: Small investment, large impact

COMPARISON OF AXILLARY AND TEMPORAL ARTERY THERMOMETRY IN PRETERM NEONATES

Abstract

The purpose of this study is to compare the accuracy of infrared temporal artery thermometry with axillary thermometry in a cohort of preterm neonates between 28 and 36 weeks postmenstrual age.

Neonates were randomly assigned to temperature measurement order (axillary followed by temporal artery or temporal artery followed by axillary). Temperature pairs were taken once during the day shift and once during the night shift. Behavioral states were assessed before, during, and after temperature measurement.

Neonates were predominantly female (64.7%) with a mean age of 6.6 days and a mean gestational age of 32.7 weeks, and most were cared for in incubators (n = 55). Noninferiority was observed between the two temperature methods (Holm-Bonferroni criterion = .025, p < .001). There was no statistically significant difference in the behavioral states of the neonates between the two temperature methods. It took nurses significantly longer to use the axillary thermometer than to use the temporal artery thermometer (p < .001).

Conclusions

- Temporal artery temperature measurements were as accurate as axillary temperature measurements in low-birth-weight neonates in the NICU.
- Nurses spent less time measuring with the temporal artery method than with the axillary method.

Title

Comparison of Axillary and Temporal Artery Thermometry in Preterm Neonates

Publication

Journal of Obstetric, Gynecologic, & Neonatal Nursing

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• /

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Preventing unplanned perioperative hypothermia in children

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AORN Journal

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Beedle, S. E., Phillips, A., Wiggins, S., & Struwe, L.

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Link

<https://aornjournal.onlinelibrary.wiley.com/doi/abs/10.1016/j.aorn.2016.12.002>

**MINIMISING THE COSTS OF TEMPERATURE MONITORING
IN HOSPITALS**

Abstract

WILL BE ADDED WHEN PDF IS READY

Conclusions

• /

Title

Minimising the costs of temperature monitoring in hospitals

Publication

Postgraduate Medical Journal

Authors

Kuymana, C.

Publication date

/

Link

<https://pmj.bmj.com/content/93/1104/580>

**TOOLKIT FOR IMPLEMENTATION OF TEMPORAL ARTERY
THERMOMETERS FOR NEONATES**

??

COST MINIMISATION ANALYSIS OF THERMOMETRY IN TWO DIFFERENT HOSPITAL SYSTEMS

Abstract

Temperature monitoring can be accomplished by various methods, including oral (OT), rectal (RT), axillary (AT), tympanic membrane (TMT) and temporal artery (TAT) thermometry, with varying amounts of cost incurred by healthcare systems.

At UHCZ, use of TAT would bring budget savings regardless of the number of devices per bed and the number of years observed. Savings would range from US\$0.08 million (one device per bed, impact for 1 year) to US\$1.8 million (one device per 10 beds, impact for 5 years). At UMH, use of TAT would lead to budget savings if one device per 10 beds were acquired, but only over a period of 3 or 5 years. Other TAT scenarios were associated with budget costs at UMH even after a period of 5 years.

Conclusion

- Sensitivity analyses showed that the price of current consumables had the highest impact on the model in both hospital settings, with TAT up to 10 times cheaper in that regard over TMT at UHCZ, potentially leading to considerable budget savings within a year of hospital-wide implementation.

Title

Cost minimisation analysis of thermometry in two different hospital systems

Publication

Postgrad Medical Journal

Authors

Hayes, K., Shepard, A., Cesarec, A., et al.

Publication date

January 2017

Link

<https://pmj.bmj.com/content/93/1104/603>

SUMMARY SAFETY REVIEW - EAR AND FOREHEAD (CONTACT) INFRARED THERMOMETERS (VARIOUS BRANDS) - ASSESSING THE POTENTIAL RISK OF INACCURACY IN CHILDREN UNDER 2 YEARS OLD

????

**PRECISION, SENSITIVITY AND PATIENT PREFERENCE OF NON-
INVASIVE THERMOMETERS IN A PEDIATRIC SURGICAL ACUTE
CARE SETTING****Abstract**

The purpose of the project was to compare the temporal artery thermometer (TAT) to the digital probe thermometer readings at axillary or oral sites, to determine the relative precision and sensitivity of the three methods of thermometry, to compare their readings to core temperature when feasible, and to survey patient and family thermometer preferences.

The study had a randomized crossover design in a 70-bed surgical unit over eight months. Two sets of temperature measurements were obtained for each patient: TAT, axillary, oral (depending on patient ability) and a bladder temperature representing core body temperature (when available). Each method was used twice on each patient, to examine within-method precision. Following measurement, patients or caregivers provided their thermometer preference. For younger/nonverbal patients, a professional observer recorded a disruption score. N=298 patients were enrolled.

TAT was more precise than oral and axillary thermometers ($p < 0.001$ vs. axillary, $p = 0.001$ vs. oral). TAT measurements were higher on average than axillary and oral, by 0.7°C and 0.6°C respectively ($p < 0.001$). TAT's disruption score for younger patients was 0.6 points lower on average than axillary ($p < 0.001$). 84% of patients and families who indicated a clear thermometry preference chose TAT. Only 3 patients had bladder-temperature devices, and therefore accuracy could not be analyzed.

Conclusions

- TAT is more precise, more fever sensitive, less disruptive to younger children, and more preferred by patients and families.
- TAT is an acceptable temperature measure that could be substituted for oral or axillary temperature in acute care pediatric settings.

Title

Precision, Sensitivity and Patient Preference of Non-Invasive Thermometers in a Pediatric Surgical Acute Care Setting.

Publication

Journal of Pediatric Nursing

Authors

Opersteny, E., Anderson, H., Bates, J., Davenport, K., Husby, J., Myking, K., & Oron, A. P.

Publication date

July/August 2017

Link

<https://www.ncbi.nlm.nih.gov/pubmed/28728767>

**TYMPANIC TEMPERATURE VERSUS TEMPORAL
TEMPERATURE IN PATIENTS WITH PYREXIA AND CHILLS****Abstract**

Accurate body temperature (BT) measurement is critical for immediate and correct estimation of core BT; measurement of changes in BT can provide physicians the initial information for selecting appropriate diagnostic approach and may prevent unnecessary diagnostic investigation. This study aimed to assess differences in tympanic and temporal temperatures among patients with fever in different conditions, especially in those with and without chills. This prospective study included patients from the emergency department between 2011 and 2012. All temperature measurements were obtained using tympanic thermometers and infrared skin thermometers. Differences in tympanic and temporal temperatures were analyzed according to 6 age groups, 5 ambient temperature groups, and 6 tympanic and temporal temperature subgroups. General linear model analysis and receiver operating characteristic curve analysis were used to estimate the differences in mean tympanic and temporal temperatures. Of the 710 patients enrolled, 246 had tympanic temperature more than 38.0°C , including 46 with chills (18.7%). Fourteen patients (3.0%) had chills and tympanic temperature less than 38°C . In the tympanic temperature subgroup of 39.0 to less than 39.5°C , approximately one-third of the patients had chills (32.3%). In the tympanic temperature subgroup of 38.0 to less than 39.0°C , the tympanic temperature was 0.4°C higher than the temporal temperature in patients without chills and 0.9°C higher in patients with chills. In the tympanic temperature subgroup of 39.0°C or more, tympanic temperature was 0.7°C higher than temporal temperature in patients without chills and 0.8°C higher in patients with chills. Temporal thermometer is more reliable in the age group of less than 1 year and 18 to less than 65 years. When the patients show tympanic temperature range of 38.0 to less than 39.0°C , 0.4°C should be added for patients without chills and 0.9°C for patients with chills to obtain core temperature. However, in patients with tympanic temperature of 39.0°C or more, 0.7°C to 0.8°C should be added, regardless of the presence of chills.

Conclusions

- [WILL BE ADDED WHEN PDF IS AVAILABLE]

Title

Tympanic temperature versus temporal temperature in patients with pyrexia and chills

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The impact of user technique on temporal artery thermometer measurements

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Nursing Critical Care

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Barry, L., Branco, J., Kargbo, N., Venuto, C., Werfel, E., Barto, D., & Glasofer, A.

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https://www.nursingcenter.com/journalarticle?Article_ID=3647718&Journal_ID=606913&Issue_ID=3647688

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Human temperatures for syndromic surveillance in the emergency department: data from the autumn wave of the 2009 swine flu (h1n1) pandemic and a seasonal influenza outbreak

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Is skin temperature measurement in PACU an accurate reflection of core temperature?

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IMPROVING PEDIATRIC TEMPERATURE MEASUREMENT IN THE ED

Abstract

ED care providers have long debated which of the various methods of temperature measurement of pediatric patients is best. While the efficacy and accuracy of temporal artery, tympanic membrane, axillary, and infrared temperature measurement have been studied, the gold standard has been rectal temperature measurement. But despite its accuracy, this method causes children with noninfectious complaints and their families unnecessary distress and adds significant time to the triage process. In response, a group of ED staff nurses at a multihospital health system conducted an evidence-based quality improvement project to determine the best practice for accurate temperature measurement in children younger than five years who presented to the ED. The project included an exhaustive literature search, a review of relevant studies, the development of a table of evidence, a presentation of the findings, and recommendations for practice change. This article describes the project and the adoption of temporal artery thermometry, a painless, noninvasive screening method that provided consistently accurate temperature measurement as well as increased patient and nurse satisfaction and a shorter triage process.

Conclusion

- WILL BE ADDED WHEN PDF IS READY

Title

Improving pediatric temperature measurement in the ed

Publication

The American Journal of Nursing

Authors

Hurwitz, B., Brown, J., & Altmiller, G.

Publication date

September 2015

Link

<https://www.ncbi.nlm.nih.gov/pubmed/26312806>

EQUIVALENCE OF TEMPERATURE MEASUREMENT METHODS IN THE ADULT HEMATOLOGY/ONCOLOGY POPULATION

Abstract

Neutropenic patients often become febrile and require frequent temperature monitoring. Monitoring core body temperature is considered the gold standard for accuracy. Taking oral temperature is the recommended noninvasive practice when core body temperature cannot be obtained; however, neutropenic patients often, for various reasons, are unable to tolerate an oral probe.

The purpose of this article is to determine the equivalence of temperatures taken via temporal artery, axillary, and oral methods, and to determine the best alternative to the oral method in the adult hematology/oncology population.

Conclusions

- The temporal method is a potential noninvasive alternative to the oral method for the adult hematology/oncology population.

Title

Equivalence of temperature measurement methods in the adult hematology/oncology population

Publication

Clinical Journal of Oncology Nursing

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Comparison of temporal to pulmonary artery temperature in febrile patients.

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Dimensions of Critical Care Nursing

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Furlong, D., Carroll, D., Finn, C., Gay, D., Gryglik, C., & Donahue, V.

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TEMPORAL ARTERY THERMOMETRY TO DETECT PEDIATRIC FEVER

Abstract

This research investigated effectiveness of temporal artery thermometry (TAT) to detect high rectal fever in children ≥ 91 days and ≤ 4 years old. Rectal temperature was initially evaluated immediately followed by TAT. As expected, the difference between mean rectal ($38.05 \pm .99$ °C) and mean TA ($37.55 \pm .8$ °C) temperatures in subjects (N = 239) was significant ($p < .0001$). Linear regression revealed TAT underestimated rectal thermometry with greater frequency at higher temperatures.

This observation provides probable explanation for the disparity between these thermometry methods. A TAT sensitivity of 75% and specificity of 85% were determined for detecting high fever (39 °C)-a finding clinically unacceptable. In contrast, among the small number of injured subjects enrolled, TAT detected high rectal fever with 100% sensitivity and specificity.

This finding, if confirmed, suggests TAT screening for well and injured children has potential for clinical practice by diminishing rectal measurements and their associated risks in the acute care and/or ambulatory practice setting.

Conclusions

- WILL BE ADDED WHEN PDF IS READY

Title

Temporal Artery Thermometry to Detect Pediatric Fever

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Tympanic, Infrared Skin, and Temporal Artery Scan Thermometers Compared with Rectal Measurement in Children: A Real-Life Assessment

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Current Therapeutic Research

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Allegaert, K., Casteels, K., van Gorp, I., & Bogaert, G.

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<https://www.ncbi.nlm.nih.gov/pmc/articles/PMC4008772/>

COMPARISON OF TEMPORAL ARTERY TEMPERATURE MEASUREMENT WITH STANDARD TEMPERATURE MEASUREMENT IN CRITICALLY ILL CHILDREN

Abstract

The objective of this study was to compare temporal artery (TA) temperature to bladder, rectal and axillary temperatures.

Sixty-nine participants aged 0 to 12 had temperature measurements from three sources: TA, bladder or rectal, and axillary. Participants were stratified based on age (< 1 year and 1–12 years) and method of temperature monitoring. Comparison of TA with other measures were computed based on deviation of $\geq 0.5^{\circ}\text{C}$.

The difference between bladder and TA temperature in <1 year group ($p < .001$) was significant at 0.55°C with bladder lower than TA temperature (95% CI: 0.24 to 0.85); however, bladder temperatures in the 1–12 year group ($p=0.57$) showed a difference of 0.1°C lower (95% CI: -0.23 to 0.42) to TA temperatures. Rectal and TA temperature did not differ significantly in either group (<1 year mean (95% CI) $=0.31^{\circ}\text{C}$ (-0.15, 0.77) ($p = .19$) and 1–12 years mean (95% CI) $=0.29^{\circ}\text{C}$ (-0.17, 0.75) ($p=0.22$)). Axillary to TA temperatures, in both age groups, reflected a significant difference (<1 year mean (95% CI) $=0.91^{\circ}\text{C}$ (0.6, 1.2) ($p < 0.001$) and 1–12 years mean (95% CI) $=0.87^{\circ}\text{C}$ (0.6, 1.2) ($p < 0.001$)).

Conclusions

- The results of this study support use of TA thermometry in the pediatric ICU as an acceptable alternative to rectal temperatures and confirm that TA thermometry is a more accurate estimate of core temperature than axillary.
- Varied bias between TA and bladder thermometry depended on the age of the patient.

Title

Comparison of temporal artery temperature measurement to standard temperature measurement in critically ill children

Publication

Pediatric Critical Care Medicine

Authors

Merrill, K., Ridling, D., Hawk, H., Jones, S., Nelson, J., Ruddy, M., Wilson, A., & Yalon, L.

Publication date

May 2014

Link

https://journals.lww.com/pccmjournal/Fulltext/2014/05001/ABSTRACT_179___COMPARISON_OF_TEMPORAL_ARTERY.176.aspx

RELIABILITY OF TEMPORAL ARTERY THERMOMETERS

Abstract

WILL BE ADDED WHEN PDF IS READY

Conclusions

- /

Title

Reliability of temporal artery thermometers

Publication

The Journal of Adult Health

Authors

McConnell, E., Senseney, D., George, S. S., & Whipple, D.

Publication date

November/December 2013

Link

<https://www.ncbi.nlm.nih.gov/pubmed/24600936>

LARYNGEAL VIEW AND TEMPERATURE MEASUREMENTS WHILE USING THE PERILARYNGEAL AIRWAY (COBRA-PLUS) IN CHILDREN

Abstract

The Cobra-PLUS perilaryngeal airway (CP) is a modification of the Cobra perilaryngeal airway. It has a distal curve for easier placement and a thermistor on the pharyngeal cuff. We assessed the orientation of the larynx to the CP and compared temperatures measured using CP to temporal arterial (TA) and infrared tympanic (T) thermometers.

American Society of Anesthesiologists (ASA) physical status 1 and 2 children 0-18 years old undergoing general anesthesia using CP were grouped into different weight cohorts. A fiberoptic scope was inserted through the CP, and laryngeal views were recorded and graded off line. Temperatures were measured from the CP, TA, and T at 15-min intervals for four readings or until the end of surgery. The CP was removed, while the patient was deeply anesthetized. Eighty subjects were analyzed. 87.5% (cohort range 75-95%) had an unobstructed view of the larynx. No serious adverse effects noted. Three hundred and sixteen temperature data points were recorded for each measured site. CP temperatures were consistently lower than TA and T with a bias of 0.9 and 0.6°C, respectively. Using temperatures measured at time 0 and 15 min, CP was associated with a larger intraclass correlation coefficient and smaller repeatability coefficient when compared to TA or T (ICC 0.65, 0.46, 0.44 and RC 0.78, 1, 1.36, respectively), indicating it had a better measure and remeasure reliability.

Conclusions

- The CP has a better orientation to the larynx compared with its previous version.
- It may be used to reliably trend intraoperative temperatures.

Title

Laryngeal view and temperature measurements while using the perilaryngeal airway (Cobra-PLUS) in children

Publication

Pediatric Anesthesia

Authors

Tan, G. M., Galinkin, J. L., Pan, Z., & Polaner, D. M.

Publication date

December 2013

Link

[Go to the full Study in PDF >](#)

STUDY RECOMMENDS USE OF PROFESSIONAL TEMPORAL THERMOMETER IN ADULTS

Abstract

Article in Oncology Nurse Advisor that states: "Using the Exergen professional TAT-5000 temporal thermometer as a noninvasive alternative to measuring temperature orally in adult hematology/oncology patients to detect febrile neutropenia is supported by research conducted at Moffitt Cancer Center, a study presented at the Oncology Nursing Society (ONS) 38th Annual Congress outlined."

Title

Study recommends use of professional temporal thermometer in adults

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标题
Comparison of rectal, axillary, tympanic, and temporal artery thermometry in the pediatric emergency room

发表期刊
Pediatric Emergency Care

作者
Batra, P., & Goyal, S.

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Link
<https://www.ncbi.nlm.nih.gov/pubmed/23283266>

COMPARISON OF TEMPORAL ARTERY TO MERCURY AND DIGITAL TEMPERATURE MEASUREMENT IN PEDIATRICS

Abstract
The aim of the study was to compare the temporal artery thermometer measurements with the mercury and digital axillary thermometer measurements in children.

This study was conducted at the Pediatric Emergency Department of Akdeniz University Hospital over a three-month period in Turkey. The sample for the study comprised 218 children (aged 0-18years). Three different methods were applied to each patient at the same time. After acquiring necessary institutional permission to conduct the study, the informed consent to participate was obtained from parents before the subjects were included in the study. The data were evaluated using general linear models. The differences between the groups were analyzed by Least Significant Difference method.

The average temperature measured by temporal artery, mercury and digital thermometers were 38.9°C, 38.3°C and 38°C respectively.

Conclusions

- Temporal artery thermometer values might be considered as core temperature.
- Rectal temperature is about 2°F (1°C) higher than an axillary temperature.
- In our study the difference between the temperature measured values was found to be consistent with the range provided in the literature.
- Temporal artery thermometers are recommended especially pediatric emergency settings, where an accurate, quick and safe body temperature measurement is of vital importance.

Title
Comparison of temporal artery to mercury and digital temperature measurement in pediatrics.

Publication
International Emergency Nursing

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Publication date
July 2014

Link
<https://www.ncbi.nlm.nih.gov/pubmed/24183491>

**ARE TEMPORAL ARTERY TEMPERATURES ACCURATE ENOUGH
TO REPLACE RECTAL TEMPERATURE MEASUREMENT IN
PEDIATRIC ED PATIENTS?**

Abstract

This study examined the accuracy of temporal artery and axillary temperatures compared with rectal temperatures in pediatric ED patients younger than 4 years.

A method-comparison study design was used to examine the agreement between a temporal artery or axillary thermometer and a nondisposable, rectal electronic thermometer, which is the clinical reference standard for temperature measurement in children. Temperatures were taken with each device in a convenience sample of stable, pediatric ED patients who were younger than 4 years. Bias and precision were calculated to quantify the differences between the 2 devices, as well as the percentage of temporal artery and axillary temperatures that were $>\pm 1.0^{\circ}\text{C}$ and $>\pm 1.5^{\circ}\text{C}$ higher or lower than the rectal temperature.

Conclusions

- A total of 52 pediatric ED patients were studied over a 10-month period.
- Bias and precision for the temporal artery and axillary devices were $-0.46^{\circ}\text{C} \pm 0.50^{\circ}\text{C}$ and $-0.93^{\circ}\text{C} \pm 0.49^{\circ}\text{C}$, respectively.
- The percentage of temporal artery and axillary temperatures that were $>\pm 1.0^{\circ}\text{C}$ and/or $>\pm 1.5^{\circ}\text{C}$ above or below the clinical reference temperature were 15% and 6%, respectively, for the temporal artery thermometer and 39% and 14%, respectively, for the axillary thermometer.

Title

Are temporal artery temperatures accurate enough to replace rectal temperature measurement in pediatric ED patients?

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/

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· /

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Temporal Artery Thermometry Use in Pediatric Patients in the Post-Anesthesia Care Unit

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ASPAN 31st National Conference 2012

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Fratto, L., Hogan, K., & Kenney, K.

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[https://www.jopan.org/article/S1089-9472\(12\)00124-4/fulltext](https://www.jopan.org/article/S1089-9472(12)00124-4/fulltext)

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Link
<https://www.ncbi.nlm.nih.gov/pmc/articles/PMC3440892/>

[ONLY THE PRINT FRIENDLY PDF OF THE ONLINE ARTIKEL AVAILABLE]

UNDERSTANDING THE DIFFERENT METHODS FOR TAKING A TEMPERATURE

Abstract
[Will be added when PDF is bought]

Conclusions

Title
Understanding the Different Methods for Taking a Temperature

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Blanket warmers revisited and temporal thermometers

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Biomedical Instrumentation and Technology

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Temporal artery thermometer use in the prehospital setting

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Canadian Journal of Emergency Medicine

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Carleton, E., Fry, B., Mulligan, A., Bell, A., & Brossart, C.

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https://www.cambridge.org/core/services/aop-cambridge-core/content/view/40F902A46E410120888FF6B183A2FF4E/S1481803500000877a.pdf/temporal_artery_thermometer_use_in_the_prehospital_setting.pdf

COMPARISON OF TEMPORAL ARTERY AND AXILLARY TEMPERATURES IN HEALTHY NEWBORNS

Abstract

The objective of this paper is to compare temperature readings of temporal artery and axillary thermometers in healthy late preterm and term infants in an effort to standardize practice.

Temperatures were taken at regular intervals per unit protocol. At each interval temporal and axillary temperatures were recorded.

Temporal temperatures were significantly higher (M = 36.9°C, SD = .59) than axillary temperatures (M = 36.7°C, SD = .68), $t(124) = 6.74$, $p < .0001$. Although statistical significance was shown between the two groups, no meaningful clinical difference was detected.

Conclusion

- Our study findings supported a new nursing practice standard for measuring infant temperatures in our mother/baby unit. Using temporal artery thermometers is now our unit's standard of care for healthy newborns.

Title

Comparison of Temporal Artery and Axillary Temperatures in Healthy Newborns

Publication

Journal of Obstetric, Gynecologic, & Neonatal Nursing

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- [WILL BE ADDED WHEN PDF IS READY]

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Agreement between temporal artery, oral, and axillary temperature measurements in the perioperative period

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Journal of PeriAnesthesia Nursing

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Barringer, L. B., Evans, C. W., Ingram, L. L., Tisdale, P. P., Watson, S. P., & Janken, J. K.

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Link

<https://www.ncbi.nlm.nih.gov/pubmed/21641529>

**EVIDENCE-BASED PRACTICE HABITS:
PUTTING MORE SACRED COWS OUT TO PASTURE**

Abstract

For excellence in practice to be the standard for care, critical care nurses must embrace evidence-based practice as the norm. Nurses cannot knowingly continue a clinical practice despite research showing that the practice is not helpful and may even be harmful to patients.

This article is based on 2 presentations on evidence-based practice from the American Association for Critical-Care Nurses' 2009 and 2010 National Teaching Institute and addresses 7 practice issues that were selected for 2 reasons. First, they are within the realm of nursing, and a change in practice could improve patient care immediately. Second, these are areas in which the tradition and the evidence do not agree and practice continues to follow tradition. The topics to be addressed are (1) Trendelenburg positioning for hypotension, (2) use of rectal tubes to manage fecal incontinence, (3) gastric residual volume and aspiration risk, (4) restricted visiting policies, (5) nursing interventions to reduce urinary catheter-associated infections, (6) use of cell phones in critical care areas, and (7) accuracy of assessment of body temperature. The related beliefs, current evidence, and recommendations for practice related to each topic are outlined.

Title

Evidence-Based Practice Habits: Putting More Sacred Cows Out to Pasture

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Clinical Nurse Specialists Empowering Staff to Improve Patient Outcomes in Temperature Measurement: From PI/EBP to Nursing Research

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National Association of Clinical Nurse Specialists National Conference Abstracts

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Brenda, A., Karen, S., & Rod, D.

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https://www.nursingcenter.com/journalarticle?Article_ID=1123716&Journal_ID=54033&Issue_ID=1123586

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**ACCURACY OF TEMPORAL ARTERY THERMOMETRY IN
NEONATAL INTENSIVE CARE INFANTS**

Abstract

The objective of this study is to determine the accuracy of temporal artery and axillary temperatures and the discomfort level of stable neonates during temperature measurement.

Bias and precision for the temporal artery and axillary devices were 0.30 ± 0.44 and 0.28 ± 0.33 , respectively. Analysis of variance found significant differences between both temporal and axillary temperatures compared to rectal temperatures ($P < .01$). Statistical differences were small and did not represent a clinically important difference. No statistical difference was found between temporal artery and axillary temperatures ($P = .81$). Increases in neonate discomfort after temperature measurement were significantly greater with axillary than increases after temporal artery temperature measurement ($P = .03$).

Conclusions

- This study found that body temperature measured with the temporal artery thermometer was similar to temperatures obtained with an axillary thermometer in stable, afebrile neonates.
- The use of temporal artery thermometry appears to be an acceptable approach for noninvasive temperature measurement in neonates, which causes less discomfort in neonates.

Title

Accuracy of temporal artery thermometry in neonatal intensive care infants.

Publication

Advances in Neonatal Care

Authors

Lee, G., Flannery-Bergey, D., Randall-Rollins, K., Curry, D., Rowe, S., Teague, M., Tuininga, C., & Schroeder, S.

Publication date

February 2011

Link

<https://www.ncbi.nlm.nih.gov/pubmed/21285659>

ACCURACY OF TEMPORAL ARTERY THERMOMETRY IN NEONATAL INTENSIVE CARE UNIT INFANTS

Abstract

The purpose of this study is to determine the accuracy of temporal artery and axillary temperatures and the discomfort level of stable neonates during temperature measurement.

Convenience sample of neonates between the ages of 32 and 40 weeks' gestation cared for in an isolette or crib.

A method-comparison design was used to compare different methods for noninvasive temperature monitoring (infra-red temporal artery; axillary electronic) to core body temperatures (indwelling rectal probe).

Bias and precision for the temporal artery and axillary devices were 0.30 ± 0.44 and 0.28 ± 0.33 , respectively. Analysis of variance found significant differences between both temporal and axillary temperatures compared to rectal temperatures ($P < .01$). Statistical differences were small and did not represent a clinically important difference. No statistical difference was found between temporal artery and axillary temperatures ($P = .81$). Increases in neonate discomfort after temperature measurement were significantly greater with axillary than increases after temporal artery temperature measurement ($P = .03$).

Conclusions

- This study found that body temperature measured with the temporal artery thermometer was similar to temperatures obtained with an axillary thermometer in stable, afebrile neonates.
- The use of temporal artery thermometry appears to be an acceptable approach for noninvasive temperature measurement in neonates, which causes less discomfort in neonates.

Title

Accuracy of temporal artery thermometry in neonatal intensive care infants.

Publication

Advances in Neonatal Care

Authors

Lee, G., Flannery-Bergey, D., Randall-Rollins, K., Curry, D., Rowe, S., Teague, M., Tuininga, C., & Schroeder, S.

Publication date

February 2011

Link

<https://www.ncbi.nlm.nih.gov/pubmed/21285659>

EFFICACY OF STANDARD DOSES OF IBUPROFEN ALONE, ALTERNATING, AND COMBINED WITH ACETAMINOPHEN FOR THE TREATMENT OF FEBRILE CHILDREN

Abstract

Many pediatricians recommend, and many parents administer, alternating or combined doses of ibuprofen and acetaminophen for fever. Limited data support this practice with standard US doses. This study compared the antipyretic effect of 3 different treatment regimens in children, using either ibuprofen alone, ibuprofen combined with acetaminophen, or ibuprofen followed by acetaminophen over a single 6-hour observation period.

Febrile episodes from children aged 6 to 84 months were randomized into the 3 treatment groups: a single dose of ibuprofen at the beginning of the observation period; a single dose of ibuprofen plus a single dose of acetaminophen at the beginning of the observation period; or ibuprofen followed by acetaminophen 3 hours later. Ibuprofen was administered at 10 mg/kg; acetaminophen at 15 mg/kg. Temperatures were measured hourly for 6 hours using a temporal artery thermometer. The primary outcome was temperature difference between treatment groups. Adverse-event data were not collected in this single treatment period study.

Sixty febrile episodes in 46 children were assessed. The mean (SD) age of the children was 3.4 (2.2) years, and 31 (51.7%) were girls. Differences among temperature curves were significant ($P < 0.001$; the combined and alternating arms had significantly better antipyresis compared with the ibuprofen-alone group at hours 4 to 6 (hour 4, $P < 0.005$; hours 5 and 6, $P < 0.001$). All but one of the children in the combined and alternating groups were afebrile at hours 4, 5, and 6. In contrast, for those receiving ibuprofen alone, 30%, 40%, and 50% had temperatures $>38.0^{\circ}\text{C}$ at hours 4, 5, and 6, respectively (hour 4, $P = 0.002$; hours 5 and 6, $P < 0.001$).

Conclusions

- During a single 6-hour observation period for these participating children, combined and alternating doses of ibuprofen and acetaminophen provided greater antipyresis than ibuprofen alone at 4 to 6 hours.

Title

Efficacy of standard doses of Ibuprofen alone, alternating, and combined with acetaminophen for the treatment of febrile children

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· /

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Body temperature surveillance and reporting in the emergency department: a practical sentinel for pandemics and bioterrorism

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Society for Academic Emergency Medicine Annual Meeting 2010

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· /

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Temperature Measurement in Patients Undergoing Colorectal Surgery and Gynecology Surgery: A Comparison of Esophageal Core, Temporal Artery, and Oral Methods

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Journal of PeriAnesthesia Nursing

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Emily, M., Sendelbach, S., Hodges, J. S., Gustafson, C., Machemer, C., Johnson, D., & Reiland, L.

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[https://www.jopan.org/article/S1089-9472\(10\)00040-7/abstract](https://www.jopan.org/article/S1089-9472(10)00040-7/abstract)

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标题
Comparison of the accuracy of body temperature measurements with temporal artery thermometer and axillary mercury thermometer in term newborns

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标题
Comparison of Temporal Artery to Rectal Temperature Measurements in Children Up to 24 Months

发表期刊
Journal of Pediatric Nursing

作者
Carr, E. A., Wilmoth, M. L., Eliades, A. B., Baker, P. J., Heisroth, K. L., & Stoner, K. H.

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THE NEED FOR AN ACCURATE NONINVASIVE THERMOMETER

Abstract

WILL BE ADDED WHEN PDF IS READY

Conclusion

- /

Title

The need for an accurate noninvasive thermometer

Publication

Anesthesia & Analgesia

Authors

Harper, C. M.

Publication date

July 2009

Link

https://www.researchgate.net/publication/26300715_The_Need_for_an_Accurate_Noninvasive_Thermometer

NONINVASIVE TEMPERATURE MONITORING IN POSTANESTHESIA CARE UNITS

Abstract

Initial postoperative core temperature is a physician and hospital performance measure. However, the extent to which core temperature changes during emergence from anesthesia and transport from the operating room to the postanesthesia care unit (PACU) remains unknown. Similarly, the accuracy of many noninvasive temperature-monitoring methods used in the PACU has yet to be quantified. This study, therefore, quantified the change in core temperature occurring during emergence and transport and evaluated the accuracy and precision of eight noninvasive thermometers in the PACU.

In 50 patients having laparoscopic surgery, the authors measured temperatures upon PACU arrival and 30 and 60 min thereafter. Monitoring methods included oral, axillary, temporal artery, forehead skin-surface, forehead liquid-crystal display, infrared aural canal, deep forehead, and deep chest. Bladder temperature was used as the reference and was also measured at the end of surgery. The primary outcome was agreement between individual temperatures from each method and bladder temperature in the PACU. A priori, the authors chose 0.5 degrees C as a clinically important temperature deviation. Bladder temperature increased 0.2 +/- 0.3 degrees C (95% confidence interval 0.1 to 0.3 degrees C), $P < 0.001$, during transport. None of the tested noninvasive thermometers was consistently within 0.5 degrees C of bladder temperature. However, oral, deep forehead, and temporal artery temperatures were significantly better than other methods and agreed reasonably well with bladder temperature.

Conclusion

- Invasive temperature monitoring available intraoperatively is more accurate than any generally available postoperative methods.
- Physician performance measures should therefore not be based exclusively on postoperative temperatures.
- Among the generally available postoperative monitoring methods, electronic oral thermometry appears to be the best.

Title

Noninvasive temperature monitoring in postanesthesia care units.

Publication

Anesthesiology

Authors

Langham, G. E., Maheshwari, A., Contrera, K., You, J., Mascha, E., & Sessler, D. I.

Publication date

July 2009

Link

<https://www.ncbi.nlm.nih.gov/pubmed/19512860>

DO OUR PATIENTS HAVE HYPOTHERMIA? TEMPORAL VERSUS ORAL THERMOMETERS

Abstract

PACU and Phase II Recovery used two different type thermometers for measurement. A great variance was noted in the results. A random sampling of 44 outpatients was used comparing the results of temporal thermometers at discharge from PACU to results of oral thermometers on admission to Phase II. The study was over a month's time using a tool developed by staff. The data showed greater than 40% had more than 1 degree variance and 11% had more than 2 degree variance.

Conclusions

[Will be added when PDF is bought]

Title

Do Our Patients Have Hypothermia? Temporal Versus Oral Thermometers

Publication

Journal of PeriAnesthesia Nursing

Authors

Pittman, R., & Waters, R.

Publication date

June 2009

Link

[https://www.jopan.org/article/S1089-9472\(09\)00225-1/abstract](https://www.jopan.org/article/S1089-9472(09)00225-1/abstract)

INFRA-RED THERMOMETRY: THE RELIABILITY OF TYMPANIC AND TEMPORAL ARTERY READINGS FOR PREDICTING BRAIN TEMPERATURE AFTER SEVERE TRAUMATIC BRAIN INJURY

Abstract

Temperature measurement is important during routine neurocritical care especially as differences between brain and systemic temperatures have been observed. The purpose of the study was to determine if infra-red temporal artery thermometry provides a better estimate of brain temperature than tympanic membrane temperature for patients with severe traumatic brain injury.

Brain parenchyma, tympanic membrane and temporal artery temperatures were recorded every 15-30 min for five hours during the first seven days after admission.

Twenty patients aged 17-76 years were recruited. Brain and tympanic membrane temperature differences ranged from -0.8 degrees C to 2.5 degrees C (mean 0.9 degrees C). Brain and temporal artery temperature differences ranged from -0.7 degrees C to 1.5 degrees C (mean 0.3 degrees C). Tympanic membrane temperature differed from brain temperature by an average of 0.58 degrees C more than temporal artery temperature measurements (95% CI 0.31 degrees C to 0.85 degrees C, P < 0.0001).

Conclusion

- At temperatures within the normal to febrile range, temporal artery temperature is closer to brain temperature than is tympanic membrane temperature.

Title

Infra-red thermometry: the reliability of tympanic and temporal artery readings for predicting brain temperature after severe traumatic brain injury

Publication

Critical Care

Authors

Kirk, D., Rainey, T., Vail, A., & Childs, C.

Publication date

2009

Link

<https://www.ncbi.nlm.nih.gov/pubmed/19473522>

TEMPORAL ARTERY THERMOMETRY UTILIZATION IN PEDIATRIC EMERGENCY CARE

Abstract

The objective of this study is to determine the effectiveness of temporal artery thermometry (TAT) as an alternative for temperature assessment of children 1 to 4 years of age in the pediatric emergency department.

Prospective cross-sectional study conducted at an urban children's hospital emergency department. TAT and rectal temperatures are compared in a convenience sample of children 1 to 4 years of age. Comparison of the temperatures is performed using Pearson correlation coefficient and regression analysis.

TAT and rectal temperatures are measured in 42 children 1 to 4 years of age. TAT predicts 83% of rectal temperatures. A receiver operating characteristic curve analysis shows that a cutoff of 37.7 degrees C or greater for fever in TAT is equivalent to rectal temperature greater than or equal to 38.3 degrees C with 100% sensitivity and 93.5% specificity.

Conclusions

- TAT is an effective screening tool in identifying fever in children 1 to 4 years of age.

Title

Temporal artery thermometry utilization in pediatric emergency care

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Noninvasive measurement of body temperature in critically ill patients

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American Association of Critical-Care Nurses

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Bridges, E., & Thomas, K.

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<http://ccn.aacnjournals.org/content/29/3/94.full.pdf+html>

MISGUIDED GUIDELINES ON NONINVASIVE THERMOMETRY

Abstract

[Will be added when PDF is bought]

Conclusions

Title

Misguided guidelines on noninvasive thermometry

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· /

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Tympanic membrane versus temporal artery temperatures of adult perianesthesia patients

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Journal of PeriAnesthesia Nursing

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Fetzer, S. J., & Lawrence, A.

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[https://www.jopan.org/article/S1089-9472\(08\)00197-4/abstract](https://www.jopan.org/article/S1089-9472(08)00197-4/abstract)

ACCURACY AND PRECISION OF NONINVASIVE TEMPERATURE MEASUREMENT IN ADULT INTENSIVE CARE PATIENTS

Abstract

The objective of this study is to determine accuracy and precision of oral, ear-based, temporal artery, and axillary temperature measurements compared with pulmonary artery temperature.

In 60 adults with cardiopulmonary disease and a pulmonary artery catheter, mean pulmonary artery temperature was 37.1 degrees C (SD 0.6 degrees C, range 35.3 degrees C-39.4 degrees C). Mean (SD) offset from pulmonary artery temperature (with the mean reflecting accuracy and SD reflecting precision) and confidence limits were 0.09 degrees C (0.43 degrees C) and -0.75 degrees C to 0.93 degrees C for oral measurements, -0.36 degrees C (0.56 degrees C) and -1.46 degrees C to 0.74 degrees C for ear measurements, -0.02 degrees C (0.47 degrees C) and -0.92 degrees C to 0.88 degrees C for temporal artery measurements, and 0.23 degrees C (0.44 degrees C) and -0.64 degrees C to 1.12 degrees C for axillary measurements. Percentage of pairs with differences greater than +/-0.5 degrees C was 19% for oral, 49% for ear, 20% for temporal artery, and 27% for axillary measurements. Intubation increased oral measurements compared with pulmonary artery temperatures (mean difference 0.3 degrees C, SD 0.3 degrees C, P = .001).

Conclusion

- Oral and temporal artery measurements were most accurate and precise.
- Axillary measurements underestimated pulmonary artery temperature.
- Ear measurements were least accurate and precise.
- Intubation affected the accuracy of oral measurements; diaphoresis and airflow across the face may affect temporal artery measurements.

Title

Accuracy and precision of noninvasive temperature measurement in adult intensive care patients.

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· /

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OTC device: temporal scanner TAT-2000C

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Journal of the American Pharmaceutical Association

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Canales, A. E.

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[https://www.japha.org/article/S1544-3191\(15\)31434-5/pdf](https://www.japha.org/article/S1544-3191(15)31434-5/pdf)

**THE STRUCTURE AND PROCESSES OF CARE DELIVERY
IMPACT POSTOPERATIVE NORMOTHERMIA**

Abstract

CAN'T BE FOUND ONLINE, DOES EXERGEN HAVE IT??

Conclusions

- /

Title

Publication

Authors

Publication date

Link

**PERIOPERATIVE TEMPERATURE AUDIT IN A LARGE
PEDIATRIC HOSPITAL**

Abstract

The combination of anesthetic-induced impairment of thermoregulatory control and exposure to a cool operating room environment makes most surgical patients

hypothermic. Infants and children cool more quickly than adults because their higher surface area to weight ratio favors heat loss. For the same reason, pediatric patients can be rewarmed faster than adults. The safe range for a child's core temperature was determined to be between 35.5 to 37.5oC.

The aim of this study was to determine the incidence of hypothermia after surgery in a large pediatric hospital. Hypothermia was defined as temperature < 35.5 oC at arrival in the recovery room.

Conclusions

- We found that 97.9% of pediatric patients undergoing surgery in our institution were normothermic at arrival in the recovery room
- Only 1.9% of our patients were found to be mildly hypothermic

Title

Perioperative Temperature Audit in a Large Pediatric Hospital

Publication

The Anesthesiology Annual Meeting 2007

Authors

Szmuk, P., Curry, B. P., Sheeran, P. W., Farrow-Gillespie, A. C., & Ezri, T.

Publication date

October 2007

Link

<http://www.asaabstracts.com/strands/asaabstracts/abstract.htm?year=2007&index=15&absnum=174>

TEMPERATURE MEASUREMENT IN CRITICALLY ILL ADULTS

Abstract

Despite increasing use of tympanic thermometers in critically ill patients who do not have a pulmonary artery catheter in place, variations in measurements obtained with the thermometers are still a problem.

The objective of this study was to compare the range of variability between tympanic and oral electronic thermometers.

The magnitude of error for the ThermoScan tympanic thermometer differed significantly from that of the Genius II tympanic thermometer and the SureTemp oral thermometer ($P < .001$). Application of the Bland and Altman method to frame the data on the basis of an accuracy tolerance zone of ± 0.5 degrees C indicated variability with both the oral and tympanic methods. The overall degree of variability was lower for the oral thermometer.

Conclusions

- Oral thermometers provide less variable measurements than do tympanic thermometers.
- Use of oral thermometry is recommended as the best practice method for temperature evaluation in critical care patients when measurement of core temperature via a pulmonary artery catheter is not possible.

Title

Temperature measurement in critically ill adults: a comparison of tympanic and oral methods.

Publication

American Journal of Critical Care

Authors

Giuliano, K. K., Giuliano, A. J., Scott, S. S., MacLachlan, E., Pysznik, E., Elliot, S., & Woytowicz, D.

Publication date

July 2000

Link

<https://www.ncbi.nlm.nih.gov/pubmed/10888148>

A BRIEF REPORT ON THE NORMAL RANGE OF FOREHEAD TEMPERATURE AS DETERMINED BY NONCONTACT, HANDHELD, INFRARED THERMOMETER

Abstract

[Will be added when PDF is bought]

Conclusions

Title

A brief report on the normal range of forehead temperature as determined by noncontact, handheld, infrared thermometer

Publication

American Journal of Infection Control

Authors

Ng, D. K., Chan, C., Chan, E. Y., Kwok, K., Chow, P., Lau, W., & Ho, J. C.

Publication date

May 2006

Link

[https://www.ajicjournal.org/article/S0196-6553\(05\)00005-2/abstract](https://www.ajicjournal.org/article/S0196-6553(05)00005-2/abstract)

**COMPARISON OF TEMPORAL ARTERY AND ORAL
TEMPERATURES IN THE EMERGENCY DEPARTMENT**

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[ONLY THE PRINT FRIENDLY PDF OF THE ONLINE ARTIKEL AVAILABLE]

COMPARISON OF TEMPORAL ARTERY THERMOMETER TO STANDARD TEMPERATURE MEASUREMENTS IN PEDIATRIC INTENSIVE CARE UNIT PATIENTS

Abstract

The objective of this paper is to determine the accuracy of noninvasive infrared temporal artery thermometry compared with rectal, axillary, and pulmonary artery catheter measurements in pediatric intensive care patients, and to determine whether temporal artery temperatures are affected by circulatory shock or by vasopressor use. We hypothesized that temporal artery temperatures do not differ from axillary and rectal temperatures in critically ill children, but temporal artery accuracy is decreased by shock or vasopressor use.

Mean bias was calculated between comparison pairs using each temperature method. Bland-Altman analysis demonstrated wide variability between methods. No significant differences in mean bias were seen between method pairs for all temperatures, but bias was significantly less in pulmonary artery catheter-rectal pairs compared with other method pairs. In febrile (> 38 degrees C) patients, bias in rectal-temporal artery and rectal-axillary was significantly greater than in temporal artery-axillary pairs ($p < .001$). Mean bias in pulmonary artery catheter-rectal pairs was also significantly smaller than in other pairs for all patients ($p = .008$) and febrile patients ($p = .049$). Presence of shock or vasopressor use did not significantly increase bias in any comparison pair. Sensitivity and specificity of both temporal artery and axillary for diagnosing fever were similar and improved with fever definition at temperatures > 38.5 degrees C.

Conclusion

- Temporal artery and axillary temperature measurements showed variability to rectal temperatures but had marked variability in febrile children.
- Neither was sufficiently accurate to recommend replacing rectal or other invasive methods.
- As temporal artery and axillary provide similar accuracy, temporal artery thermometers may serve as a suitable alternative for patients in whom invasive thermometry is contraindicated.

Title

Comparison of temporal artery thermometer to standard temperature measurements in pediatric intensive care unit patients

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Newborn temperature during skin-to-skin breastfeeding in couples having breastfeeding difficulties

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Birth

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Chiu, S. H., Anderson, G. C., & Burkhammer, M. D.

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<https://www.ncbi.nlm.nih.gov/pubmed/15918868>

TEMPORAL SCANNER THERMOMETRY: A NEW METHOD OF CORE TEMPERATURE MEASUREMENT IN INTENSIVE CARE PATIENTS

Abstract

Temperature measurement is a routine task of patient care, with considerable clinical impact, especially in the ICU. This study was conducted to evaluate the accuracy and variability of the Temporal Artery Thermometer (TAT) in ICU patients. Therefore, a convenience sample of 57 adult patients, with indwelling pulmonary artery catheters (PAC) in a 40-bed intensive care unit in a university teaching hospital was used.

The study design was a prospective, descriptive comparative analysis. Body temperature was thereby measured simultaneously with the TAT and the Axillary Thermometer (AT), and was compared with the temperature recording of the PAC. The use of vasoactive medication was recorded.

Mean temperature of all measurements was: PAC: 37.1 degrees C (SD: 0.87), TAT 37.0 degrees C (SD: 0.68) and axillary thermometer: 36.6 degrees C (SD: 0.94). The measurements of the TAT and the PAC were not significantly different (mean difference: 0.14 degrees C; SD: 0.51; $p = 0.33$); whereas the measurements of the PAC and the AT differed significantly (mean difference: 0.46 degrees C; SD: 0.39; $p < 0.001$). Mean difference in PAC versus TAT analyses, between patients with vasopressor therapy (0.12 degrees C; SD: 0.55), and without vasopressor therapy (0.19 degrees C; SD: 0.48) was not statistically significant ($p = 0.47$).

Conclusions

- We can conclude that the temporal scanner has a relatively good reliability with an acceptable accuracy and variability in patients with normothermia.
- The results are comparable to those of the AT, but they do not seem to be sufficient to prove any substantial benefit compared to rectal, oral or bladder thermometry.

Title

Temporal scanner thermometry: a new method of core temperature estimation in ICU patients.

Publication

Scottish Medical Journal

Authors

Myny, D., De Waele, J., Defloor, T., Blot, S., & Colardyn, F.

Publication date

February 2005

Link

<https://www.ncbi.nlm.nih.gov/pubmed/15792381>

COMPARISON OF THE TEMPORAL ARTERY AND RECTAL THERMOMETRY IN CHILDREN IN THE EMERGENCY DEPARTMENT

Abstract

Rectal thermometry, the criterion standard of temperature measurement in young children, has numerous disadvantages. This study examined the agreement between rectal versus a new temporal artery professional model (TAPM) thermometer and rectal versus a home device temporal artery consumer model (TACM) thermometer, investigated if the TAPM can safely screen for rectal fever, and determined if parents can detect rectal fever using the TACM. DESIGN, OUTCOME MEASURES, AND SUBJECTS: In this cross-sectional agreement emergency department study, 327 children <24 months of age had their temperature measured rectally and by the TAPM and TACM by a single nurse and using the TACM by the parents. Agreements were analyzed by the Bland Altman plots. Temperature cutoff to detect rectal fever $> \text{ or } = 38.0$ degrees C and $> \text{ or } = 38.3$ degrees C with sensitivities of $> \text{ or } = 90\%$ and $> \text{ or } = 95\%$, respectively, was determined for the TAPM.

The mean difference between the rectal minus TAPM was -0.19 degrees C ± 0.66 degrees C, and minus the TACM home device, it was $+0.11$ degrees C ± 0.66 degrees C. The sensitivities of TAPM temperature of $> \text{ or } = 37.7$ degrees C to detect rectal fever $> \text{ or } = 38.0$ degrees C and $> \text{ or } = 38.3$ degrees C were 90% (95% confidence interval: 0.83; 0.94) and 97% (95% confidence interval: 0.92; 0.99), respectively. The parents detected 67% and 73% of rectal fevers 38.0 degrees C and $> \text{ or } = 38.3$ degrees C, respectively.

Conclusions

- The TAPM thermometer cannot replace the rectal. However, TAPM temperature of < 37.7 degrees C can be safely used as a screen to exclude rectal fever $> \text{ or } = 38.3$ degrees C in infants 3 to 24 months of age.
- The TACM home device has insufficient ability to detect rectal fever.
- A multicenter trial is needed to validate these results across multiple emergency departments and numerous observers.

Title

Comparison of the temporal artery and rectal thermometry in children in the emergency department

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标题
When body temperature changes, does rectal temperature lag?

发表期刊
Journal of Pediatrics

作者
Greenes, D. S., & Fleisher, G. R.

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链接
[https://www.jpeds.com/article/S0022-3476\(04\)00167-2/pdf](https://www.jpeds.com/article/S0022-3476(04)00167-2/pdf)

NON-INVASIVE TEMPORAL ARTERY THERMOMETRY: PHYSICS, PHYSIOLOGY, AND CLINICAL ACCURACY, PRESENTED AT MEDICAL THERMOMETRY FOR SARS DETECTION

Abstract
Temporal artery (TA) thermometry was developed in answer to requests by pediatricians for a replacement for: 1) ear thermometry due to inaccuracy; and 2) rectal thermometry due to parents' (and most clinicians') growing dislike of the method. The underlying technology development spans some 20 years, borrowing heavily from methods invented for industrial processes and medical research. Although the forehead has been used since antiquity to detect fever, its accuracy had always been questionable until physiological artifacts were understood and overcome, and mathematical modeling of arterial heat balance at the skin has made it possible to produce accurate core temperatures entirely non-invasively with just a scan of the forehead. Clinical studies have been conclusive as to TA superiority to ear thermometry, and well on the way to being conclusive as to TA at least as accurate as rectal. The physics are relatively straightforward, but the physiological requirements are not. Underlying physiological artifacts cause errors of more than 2 deg C in non-invasive thermometry and must be reduced by an order of magnitude to provide medically useful temperatures. Patented TA technology incorporates methods of dealing with physiological artifacts to overcome these errors. Mass screening for SARS containment with this method is examined.

Conclusions

- The TA thermometer was validated by initial studies performed by qualified clinical researchers in prestigious institutions.
- Reproducibility of both the instrument and physiological assumptions was established by comparing paired left-right readings.
- Exergen and its clinical research partners concluded that the TA method was indeed suitable for routine clinical use.

Title
Non-invasive temporal artery thermometry: Physics, physiology, and clinical accuracy

Publication
Proceedings of SPIE - The International Society for Optical Engineering

Authors
Pompei, P., & Pompei, M.

Publication date
April 2004

Link
https://www.researchgate.net/publication/228597147_Non-invasive_temporal_artery_thermometry_Physics_physiology_and_clinical_accuracy

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Comparison of temporal artery, rectal and esophageal core temperatures in children: Results of a pilot study

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Journal of Pediatric and Child Health

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Al-Mukhaizeem, F., Allen, U., Komar, L., Naser, B., Roy, L., Stephens, D., Read, S., Kim, C., & Schuh, S.

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- /

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A comparison of measurements from a temporal artery thermometer and a pulmonary artery catheter thermistor

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American Journal of Critical Care

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Furlong, D., Carroll, D. L., Finn, C., Gay, D., Gryglik, C., & Donahue, V.

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<https://www.ncbi.nlm.nih.gov/pubmed/25470268>

NEW PRODUCT REVIEW: TEMPORAL ARTERY THERMOMETRY

Abstract

[WILL BE ADDED WHEN PDF IS BOUGHT]

Conclusions

Title

New product review: temporal artery thermometry

Publication

Journal of PeriAnesthesia Nursing

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Detecting fever in young infants: reliability of perceived, pacifier, and temporal artery temperatures in infants younger than 3 months of age

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Pediatric Emergency Care

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Callanan, D.

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<https://www.ncbi.nlm.nih.gov/pubmed/12972820>

TEMPORAL ARTERY TEMPERATURE MEASUREMENTS IN HEALTHY INFANTS, CHILDREN, AND ADOLESCENTS

Abstract

A noninvasive temporal artery thermometer that uses arterial heat balance technology has been compared to rectal and ear thermometry and is available in the marketplace. This study was undertaken to establish mean temperatures and temperatures 2 standard deviations above the mean for healthy infants, children, and adolescents. Temperatures were measured in healthy patients 0 to 18 years of age using a noninvasive temporal artery thermometer. Temperatures were measured in 2,346 patients. Mean temperatures and temperatures 2 standard deviations above the mean were: 37.1 degrees C (38.1 degrees C) for 383 infants 0 to 2 months; 36.9 degrees C (37.9 degrees C) for 860 children 3 to 47 months; 36.8 degrees C (37.8 degrees C) for 680 children 4 to 9 years; and 36.7 degrees C (37.8 degrees C) for 423 adolescents 10 to 18 years. There were no significant differences in temperatures in white compared to African-American children, children with or without perspiration on their forehead, or between measurements taken on the left compared to the right side of the forehead. This study provides information about temporal artery temperatures in healthy infants and children that can serve as a basis for interpreting temperature measurements in ill children when the same instrument is used.

Conclusions

- [WILL BE ADDED WHEN PDF IS BOUGHT]

Title

Temporal artery temperature measurements in healthy infants, children, and adolescents

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**COMPARISON OF THE TEMPORAL ARTERY AND RECTAL
THERMOMETRY IN CHILDREN IN THE EMERGENCY
DEPARTMENT**

[THESE LOOK SIMILAIR]

Number 65:

Schuh S, Komar L, Stephens D, Chu L, Read S, Allen U (University of Toronto/Hospital for Sick Children). Comparison of the temporal artery and rectal thermometry in children in the emergency department. Pediatric Academic Societies Annual Meeting, May 3-6, 2003, Seattle, WA.

Number 66:

Schuh S, Komar L, Stephens D, Chu L, Read S, Allen U (University of Toronto/Hospital for Sick Children). Comparison of the temporal artery and rectal thermometry in children in the emergency department. Pediatric Emergency Care, Vol 20, No. 11, Nov 2004

INSUFFICIENCY IN THERMOMETER DATA

Abstract

[Will be added when PDF is bought]

Conclusions

Title

Insufficiency in thermometer data

Publication

Anesthesia & Analgesia

Authors

Pompei, F.

Publication date

January 2009

Link

https://journals.lww.com/ccmjournals/Citation/2009/01000/Misguided_guidelines_on_noninvasive_thermometry.89.aspx#print-article-link

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标题
Infrared temporal temperature measurement

发表期刊
Journal of the Norwegian Medical Association/Tidsskriftet Den Norske Lægeforening

作者
Dybwik, K., & Nielsen, E. W.

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COMPARISON OF TEMPLE TEMPERATURES WITH RECTAL TEMPERATURES IN CHILDREN UNDER TWO YEARS OF AGE

Abstract
We assessed the agreement between rectal and noninvasive temporal artery temperature measurements in infants and children. We also evaluated the temple thermometer as a screening tool for rectal fever in this age group. Finally, we compared the performance of parents with that of nurses in using the temple thermometer. The 95% limits of agreement between the difference in rectal and average temple temperature were -1.03 and +1.52 degrees C. Mean temple temperatures obtained by parents and by nurses were similar (95% limits of agreement, -0.6 degrees C to +0.7 degrees C). A maximum temple temperature cutoff of 37.2 degrees C (99.0 degrees F) distinguished children with rectal fever of > or =38.0 degrees C with 91% sensitivity and 53% specificity. A cutoff of 37.8 degrees C (100.0 degrees F) distinguished moderate rectal fevers (> or =38.5 degrees C) with 97% sensitivity and 84% specificity. A cutoff of 38.3 degrees C (101.0 degrees F) distinguished a high rectal fever (> or =39.0 degrees C) with a sensitivity of 95% and specificity of 95%. In conclusion, temple temperatures do not reliably predict rectal temperatures, but the temple thermometer can be used as an effective screen for clinically important rectal fever in children 3-24 months old. The findings do not support use of temple temperatures to screen young infants for rectal fever > or =38.0 degrees C. Temperatures obtained by parents were comparable to those obtained by nurses.

Conclusions

- [\[WILL BE ADDED WHEN PDF IS AVAILABLE\]](#)

Title
Comparison of temple temperatures with rectal temperatures in children under two years of age

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Accuracy of a noninvasive temporal artery thermometer for use in infants

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Archives of Pediatrics & Adolescent Medicine

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Greenes, D. S., & Fleisher, G. R.

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标题
Thermometry in children

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链接
[Go to full White Paper report in PDF > \[Update link when server is ready\]](#)

**MEASURING BODY TEMPERATURE:
A SIMPLE, EFFECTIVE AND RELIABLE SOLUTION**

Abstract
....

Conclusion
•

Title
Measuring body temperature:A simple, effective and reliable solution

Publication
UMCG

Authors
??

Publication date
June 2017

Link
[Go to full White Paper report in PDF >](#)

SMALL INVESTMENT, LARGE IMPACT

Abstract

....

Conclusion

•

Title

Small investment, large impact

Publication

??

Authors

??

Publication date

??

Link

[Go to the full Case Study in PDF >](#)

**EXERGEN MEDICAL FORMS STRATEGIC DISTRIBUTION
PARTNERSHIP WITH INTERNATIONAL MEDICAL TRADE
CONSULTANTS OF ISRAEL**

Abstract

....

Conclusion

•

Title

Exergen Medical forms strategic distribution partnership with international medical trade consultants of Israel

Publication

??

Authors

??

Publication date

??

Link

[Go to the full Case Study in PDF >](#)

***AFTER UMCG, BEATRIXOORD, CENTRE FOR REHABILITATION
IN HAREN, ALSO CHOOSES INFRARED THERMOMETERS
FROM EXERGEN***

Abstract

....

Conclusion

•

Title

After UMCG, Beatrixoord, Centre for Rehabilitation in Haren, also chooses infrared thermometers from Exergen

Publication

??

Authors

??

Publication date

??

Link

[Go to the full Case Study in PDF >](#)

***UL GRANTS EXERGEN CORPORATION ISO 13485:2016
CERTIFICATE***

Abstract

....

Conclusion

•

Title

UL grants Exergen Corporation ISO 13485:2016 certificate

Publication

??

Authors

??

Publication date

??

Link

[Go to the full Case Study in PDF >](#)

**WHY SEATTLE CHILDREN'S HOSPITAL DECIDED TO ADOPT
TEMPORAL ARTERY THERMOMETERS**

Abstract

....

Conclusion

•

Title

Why Seattle Children's Hospital decided to adopt temporal artery thermometers

Publication

??

Authors

??

Publication date

??

Link

[Go to the full Case Study in PDF >](#)

**90 PERCENT OF NURSES PREFER THERMOMETER
CONNECTED TO VITAL SIGNS MONITOR**

Abstract

....

Conclusion

•

Title

90 percent of nurses prefer thermometer connected to vital signs monitor

Publication

??

Authors

??

Publication date

??

Link

[Go to the full Case Study in PDF >](#)

UNIVERSITY MEDICAL CENTRE GRONINGEN CHOOSES TAT-5000 OF EXERGEN MEDICAL AS STANDARD THERMOMETER THROUGHOUT HOSPITAL

Abstract

....

Conclusion

•

Title

University Medical Centre Groningen chooses TAT-5000 of Exergen Medical as standard thermometer throughout hospital

Publication

??

Authors

??

Publication date

??

Link

[Go to the full Case Study in PDF >](#)

TEMPORAL ARTERY THERMOMETRY VERY COMFORTABLE FOR PRETERM NEONATES

Abstract

....

Conclusion

•

Title

Temporal Artery Thermometry very comfortable for preterm neonates

Publication

??

Authors

??

Publication date

??

Link

[Go to the full Case Study in PDF >](#)

APPEALS COURT UPHOLDS EXERGEN TEMPORAL ARTERY THERMOMETER PATENT INFRINGEMENT BY HELEN OF TROY'S BRAUN AND VICKS FOREHEAD THERMOMETERS

Abstract

....

Conclusion

•

Title

Appeals Court Upholds Exergen Temporal Artery Thermometer Patent Infringement by Helen of Troy's Braun and Vicks Forehead Thermometers

Publication

??

Authors

??

Publication date

??

Link

[Go to the full Case Study in PDF >](#)

HOSPITALS SAVE 1,4 MILLION EURO BY ADOPTING TAT- 5000 TEMPORAL ARTERY THERMOMETER

Abstract

....

Conclusion

•

Title

Hospitals Save 1,4 million Euro By Adopting TAT- 5000 Temporal Artery Thermometer

Publication

??

Authors

??

Publication date

??

Link

[Go to the full Case Study in PDF >](#)