# EXERGEN

C O R P O R A T I O N

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**Competency Assessment for Exergen TA Thermometer**

Print Name \_\_\_\_\_ Date \_\_\_\_\_  
Unit \_\_\_\_\_

Choose the best answer and write the letter of that answer on the line in front of the number.

- Temporal artery thermometers measure the patients'
  - Ambient air temperature.
  - Core body temperature.
  - Skin temperature.
  - Oxid temperature.
- Core temperature will show a spike in patient temperature:
  - One to two hours later than rectal temperature.
  - At the same time as rectal temperature.
  - One to two hours sooner than rectal temperature.
- The temporal artery thermometer measures the temperature of the temporal and carotid arteries, reflecting the core temperature at the heart. In the case of the patient who has been febrile, and the fever is now breaking, the temporal artery scanner may read:
  - Lower than a rectal thermometer because the core temperature will reflect the change more rapidly than the rectum.
  - Higher than a rectal thermometer because the core temperature takes longer than the rectum to reflect the change.
  - The same because it does not matter how or where temperature is measured.
  - Lower, because environmental factors will always affect core temperature.
- The temporal artery thermometer may give inaccurate reading if:
  - The lens is dirty.
  - The side of the forehead measured has been resting on the pillow.
  - The patient has just finished drinking iced water.
  - 1 only
  - All of the above
  - 1 and 3
  - 1 and 2
  - None of the above
- Core temperature measurement reflects changes in body temperature oral or rectal temperature measurement.
  - Slower than
  - The same as
  - More quickly than
- It is important to clean the lens in the center of the probe with a cotton-tipped swab applicator (Q-Tip) dampened with an alcohol prep pad.
  - Every two weeks
  - After each use
  - When the patient is discharged
  - If lens is not shiny and mirror-like
    - 1 only
    - All of the above
    - 1 and 3
    - 1 and 4
    - None of the above

SKILL PERFORMANCE	YES	NO
1. Explains procedure to patient.		
2. Cleans probehead of temporal artery scanner between patients.		
3. Places probe flush on center of forehead, depresses button, and holds button depressed the entire time.		
4. Slides probe in a <u>straight</u> line across the forehead to the hairline.		
5. Lifts probe from forehead and touches the little soft depressor on the neck behind the earlobe.		
6. Reflexes button and reads temperature.		
7. Records temperature on bedside graphic sheet.		
8. Verbalizes intervention for any abnormal results.		
9. Demonstrates cleaning the lens of the temporal artery thermometer.		

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1. Adult - Artz, B.A., March, K.S., and Grim, R.D. Clinical nurse specialists empowering staff to improve patient outcomes in temperature measurement: From PI/EBP to nursing research. In Patton, S. P., and Moody, R. 2011 National Association of Clinical Nurse Specialists National Conference Abstracts. *Clinical Nurse Specialist: The Journal for Advanced Nursing Practice*.
2. Adult - Espenhein A (County Hospital in Herlev, Denmark). Temporal temperature measurement. *Sygeplejersken* 2006;(17):50-2.
3. Adult CCU - Bridges E, Thomas K (University of Washington). Noninvasive measurement of body temperature in critically ill patients. *Crit. Care Nurse*. 2009; 29(3): p. 94-97
4. Adult CCU - Carroll D, Finn C, Gill S, et al (Massachusetts General Hospital). A comparison of measurements from a temporal artery thermometer and a pulmonary artery catheter thermometer. *Am J Crit Care*. 2004;13:258.
5. Adult CCU - Furlong D, Carroll D, Finn C, Gay D, Gryglik C, Donahue V (2015). Comparison of Temporal to Pulmonary Artery Temperature in Febrile Patients. *Dimensions of Critical Care Nursing*. 2015 Jan-Feb; 34(1):47-52. doi: 10.1097/DCC.0000000000000090.
6. Adult CCU - Barry L, Branco J, et al. The impact of user technique on temporal artery thermometer measurements. *Nursing Critical Care: September 2016 - Volume 11 - Issue 5 - p 12 - 14*.
7. Adult CCU - Makic MB, VonRueden KT, Rauen CA, Chadwick J. Evidence-based practice habits: putting more sacred cows out to pasture. *Crit Care Nurse*. 2011 Apr;31(2):38-61; quiz 62.
8. Adult CCU - Lawson L, Bridges E, Ballou I, Eraker R, Greco S, Shively J, Sochulak V. (University of Washington). Accuracy and precision of noninvasive temperature measurement in adult intensive care patients. *Am. J. Crit. Care.*, Sep 2007; 16:5, 485-496.
9. Adult CCU - Lawson L, Bridges E, Ballou I, Eraker R, Greco S, Shively J, Sochulak V. (University of Washington). Temperature measurement in critically ill adults. *Am. J. Crit. Care.*, May 2006; 15: 324 - 346.
10. Adult ED - Bordonaro SF, McGillicuddy DC, Pompei F, Burmistrov D, Harding C, Sanchez LD. Human temperatures for syndromic surveillance in the emergency department: data from the autumn wave of the 2009 swine flu (H1N1) pandemic and a seasonal influenza outbreak. *BMC Emerg Med*. 2016 Mar 9;16(1):16. doi: 10.1186/s12873-016-0080-7.
11. Adult ED - Foy S, McGillicuddy D, Pompei F, Sanchez L (Beth Israel Medical Center, Boston MA).  
Body Temperature Surveillance and Reporting in the Emergency Department: A Practical Sentinel for Pandemics and Bioterrorism. Presented at Society for Academic Emergency Medicine Annual Meeting, Phoenix AZ , June 3-6, 2010.
12. Adult ED - Routhier D, Hostler D, Wolfson A, Wheeler M, Reynolds J (University of Pittsburgh). Comparison of temporal artery and oral temperatures in the emergency department. *ACAD EMERG MED*, May 2006, Vol. 13, No. 5, Suppl. 1, www.aemj.org , p. S99
13. Adult ICU - Dybwik K, Nielsen EW. Infrared temporal temperature measurement. *Journal of the Norwegian Medical Association* 2003; 123: 3025-6.

14. Adult ICU - Myny D, DeWaele J, Defloor T, Blot S, Colardyn F (Ghent University Hospital, Ghent, Belgium). Temporal scanner thermometry: a new method of core temperature measurement in intensive care patients. *SMJ* 2005 45(1): 15-18.
15. Adult ICU - Kirk D, Rainey T, Vail A, Childs C (University of Manchester, Salford Royal Foundation Trust). Infra-red thermometry: the reliability of tympanic and temporal artery readings for predicting brain temperature after severe traumatic brain injury. *Crit Care*. 2009 May 27;13(3):R81.
16. Adult oncology - Hughes D. Study recommends use of professional temporal thermometer in adults. *Oncology Nurse Advisor* April 27, 2013.
17. Adult oncology - Mason TM, Reich RR, et al. Equivalence of temperature measurement methods in the adult hematology/oncology population. *Clin J Oncol Nurs*. 2015 Apr;19(2):E36-40. doi: 10.1188/15.CJON.E36-E40.
18. Adult PACU - Fetzer SJ, Lawrence A (Southern New Hampshire Medical Center). Tympanic membrane versus temporal artery temperatures of adult perianesthesia patients. *J Perianesth Nurs*. 2008 Aug;23(4):230-6.
19. Adult Surgery - Calonder EM, Sendelbach S, Hodges JS, Gustafson C, Machemer C, Johnson D, Reiland L (Abbott Northwestern Hospital). Temperature measurement in patients undergoing colorectal surgery and gynecology surgery: a comparison of esophageal core, temporal artery, and oral methods. *Journal of PeriAnesthesia Nursing*, Volume 25, Issue 2, April 2010, Pages 71-78
20. Adult Surgery - McConnell E, Senseney D, George S, Whipple D. Reliability of temporal artery thermometers. *Medsurg Nursing* 2013, Nov-Dec 2013, Vol. 22/No. 6, p387
21. Ambulance - Carleton E, Fry B, Mulligan A, Bell A, Brossart C. Temporal artery thermometer use in the prehospital setting. *Canadian Journal of Emergency Medicine* 2012;14(1):7-13.
22. Biomedical - Crossley B. Blanket warmers revisited and temporal thermometers. *Biomedical Instrumentation and Technology*, March/April 2012 p147.
23. Cost Reduction Using TAT - Hayes K, Shepard A, Cesarec A, et al. Cost minimisation analysis of thermometry in two different hospital systems. *Postgrad Med J* Published Online First: 18 January 2017, doi:10.1136/postgradmedj-2016-134630
24. Cost Reduction Using TAT - Kumana C. Minimizing the costs of temperature monitoring in hospitals. *Postgrad Med J* Published Online First: 1 February 2017 doi:10.1136/postgradmedj-2017-134795
25. Methods Paper - Pompei F, Pompei M. Non-invasive temporal artery thermometry: Physics, Physiology, and Clinical Accuracy, presented at Medical Thermometry for SARS Detection, SPIE Defense and Security Symposium, available in Conference Proceedings, April, 2004.
26. Methods Paper - Pompei F. Insufficiency in thermometer data. *Anesth Analg*. 2003 Mar;96(3):908-9.
27. Methods Paper - Pompei F. RE: A brief report on the normal range of forehead temperature as determined by noncontact, handheld, infrared thermometer. *Am J Infect Control*. 2006 May;34(4):248-9.
28. Methods Paper - Pompei F. Misguided guidelines on noninvasive thermometry. *Crit Care Med*. 2009 Jan;37(1):383; author reply 383-4.
29. Neonates - Haddad, L., Smith, S., Phillips, K.D., and Heidel, R.E. (2012). Comparison of temporal artery and axillary temperatures in healthy newborns. *Journal of Obstetric, Gynecologic, & Neonatal Nursing*, 41, 383-388; doi: 10.1111/j.1552-6909.2012.01367.x
30. Neonates - Chiu SH, Anderson GC, Burkhammer MD (University of Akron/Case Western Reserve University). Newborn temperature during skin-to-skin breastfeeding in couples having breastfeeding difficulties. *Birth*. 2005 Jun;32(2):115-21.

33. NICU - Burdjalov VF, Combs A, Nachman S, Baumgart S (SUNY at Stony Brook). Non-Invasive infrared temperature assessment of the temporal artery for core temperature determination in premature neonates, Presented American Pediatric Society and the Society for Pediatric Research, May 1, 2001.
34. NICU - Lee G, Flannery-Bergey D, Randall-Rollins K, Curry D, Rowe S, Teague M, Tuininga C, Schroeder S (Exempla Lutheran Medical Center). Accuracy of temporal artery thermometry in neonatal intensive care infants. *Advances in Neonatal Care*, Vol. 11, No. 1, pp. 62-70, Feb 2011.
35. NICU - Rollins K, Flannery-Bergey D. Accuracy of temporal artery thermometry in neonatal intensive care unit infants. *JOGNN*, 40, S85-S119; 2011. DOI: 10.1111/j.1552-6909.2011.01243.x
36. NICU - Smith et al. Comparison of axillary and temporal artery thermometry in preterm neonates. *J Obstet Gynecol Neonatal Nurs*. 2018 Apr 3. pii: S0884-2175(18)30052-2. doi: 10.1016/j.jogn.2018.02.013. [Epub ahead of print]
37. PACU - Barringer LB, Evans CW, Ingram LL, Tisdale PP, Watson SP, Janken JK (Presbyterian Hospital Matthews). Agreement between temporal artery, oral, and axillary temperature measurements in the perioperative period. *J Perianesth Nurs*. 2011 Jun;26(3):143-50.
38. PACU - Bradley SL, Kwater AP, et al. Is skin temperature measurement in PACU an accurate reflection of core temperature? ASA Abstract A3182, the Anesthesiology Annual Meeting 2016, <http://www.asaabstracts.com>
39. PACU - Langham GE, Maheshwari A, Contrera K, You J, Mascha E, Sessler DI (Case Western Reserve University). Noninvasive temperature monitoring in postanesthesia care units. *Anesthesiology*. 2009 Jul;111(1):90-6. doi: 10.1097/ALN.0b013e3181a864ca
40. PACU - Martinez EA, Krenzischek D, Hobson D, Hunt D (Johns Hopkins Medical Institutions). The structure and processes of care delivery impact postoperative normothermia. *Anesthesiology* 2007; 107: A496.
41. PACU - Pittman R and Waters R (CaroMont Health Care, Gastonia, NC). Do our patients have hypothermia? Temporal versus oral thermometers. *Journal of PeriAnesthesia Nursing* Volume 24, Issue 3, June 2009, Page e18.
42. PACU - Sandlin D (Southern Hills Medical Center, Nashville TN). New Product Review: Temporal Artery Thermometry, *Journal of PeriAnesthesia Nursing*, Vol. 18, No 6 (December) 2003, pp 419-421.
43. Pediatric CCU - Operstény, Esther et al. Precision, sensitivity and patient preference of non-invasive thermometers in a pediatric surgical acute care setting. *Journal of Pediatric Nursing: Nursing Care of Children and Families* , 2017, Volume 35 , 36 - 41.
44. Pediatric CCU - Hebbar K, Fortenberry JD, Rogers K, Merritt R, Easley K. (Children's Healthcare of Atlanta at Egleston). Comparison of temporal artery thermometer to standard temperature measurements in pediatric intensive care unit patients. *Pediatr Crit Care Med*. 2005 Sep;6(5):557-61.
45. Pediatric CCU - Merrill, K. (Seattle Children's Hospital). Comparison of temporal artery temperature measurement with standard temperature measurement in critically ill children. *American Journal of Critical Care*. 2014, May, 23(3), e23.
46. Pediatric ED - Greenes DS, Fleisher GR. (Boston Childrens Hospital and Harvard Medical School). When body temperature changes, does rectal temperature lag? *Journal of Pediatrics*, 02.037, pp 824-826, September 2004.
47. Pediatric ED - Greenes DS, Fleisher GR. (Boston Childrens Hospital and Harvard Medical School). Accuracy of a noninvasive temporal artery thermometer for use in infants. *Arch Pediatr Adolesc Med*, Vol 155, pp 376-381, Mar 2001

48. Pediatric ED - Batra P, Saha A, Faridi MM. Thermometry in children. *J Emerg Trauma Shock*. 2012 Jul;5(3):246-9.
49. Pediatric ED - Hurwitz B1, Brown J, Altmiller G. Improving pediatric temperature measurement in the ED. *Am J Nurs*. 2015 Sep;115(9):48-55. doi:10.1097/01.NAJ.0000471249.69068.73.
50. Pediatric ED - Moore AH, Carrigan JD, Solomon DM, Tart RC. Temporal artery thermometry to detect pediatric fever. *Clin Nurs Res*. 2015 Oct;24(5):556-63. doi: 10.1177/1054773814557481. Epub 2014 Nov 14.
51. Pediatric ED - Reynolds M, et al. Are temporal artery temperatures accurate enough to replace rectal temperature measurement in pediatric ED patients? *J Emerg Nurs*. 2012 Nov 8. pii: S0099-1767(12)00329-7. doi: 10.1016/j.jen.2012.07.007. [Epub ahead of print]
52. Pediatric ED - Schuh S, Komar L, Stephens D, Chu L, Read S, Allen U (University of Toronto/Hospital for Sick Children). Comparison of the temporal artery and rectal thermometry in children in the emergency department. Pediatric Academic Societies Annual Meeting, May 3-6, 2003, Seattle, WA.
53. Pediatric ED - Schuh S, Komar L, Stephens D, Chu L, Read S, Allen U (University of Toronto/Hospital for Sick Children). Comparison of the temporal artery and rectal thermometry in children in the emergency department. *Pediatric Emergency Care*, Vol 20, No. 11, Nov 2004
54. Pediatric ED - Szmuk P, Curry BP, Sheeran PW, Farrow-Gillespie AC, Ezri T (UT Southwestern and Children's Medical Center, Dallas, Texas). Perioperative temperature audit in a large pediatric hospital. *Anesthesiology* 2007; 107: A1612.
55. Pediatric ED - Titus MO, Hulsey T, Heckman J, Losek JD (Medical University of South Carolina and Children's Hospital). Temporal artery thermometry utilization in pediatric emergency care. *Clinical Pediatrics*, Mar 2009; vol. 48: pp. 190 - 193.
56. Pediatric ED - Isler A., et al. Comparison of temporal artery to mercury and digital temperature measurement in pediatrics. *Int. Emerg.Nurs.* (2013) <http://dx.doi.org/10.1016/j.ienj.2013.09.003>
57. Pediatric ED - Batra P, Goyal S. Comparison of rectal, axillary, tympanic, and temporal artery thermometry in the pediatric emergency room. *Pediatr Emerg Care*. 2013 Jan;29(1):63-6. doi: 10.1097/PEC.0b013e31827b5427.
58. Pediatric PACU - Beedle SE, Phillips A, et al. Preventing unplanned perioperative hypothermia in children. *AORN J*. 2017 Feb;105(2):170-183. doi: 10.1016/j.aorn.2016.12.002.
59. Pediatric PACU - Fratto L, Hogan K, Kenney K. Temporal artery thermometry use in pediatric patients in the post-anesthesia care unit. 2012 Research and EBP Abstracts ASPAN's 31st National Conference April 15-19, 2012, Orlando, FL Inf...2012 Research and EBP Abstracts - 8/1/2012 12:43:08 PM.
60. Pediatric Surgery - Tan GM, Galinkin JL, Pan Z, Polaner DM. Laryngeal view and temperature measurements while using the perilaryngeal airway (Cobra-PLUSTM) in children. *Pediatric Anesthesia* 2013, Dec; 23(12):1180-6. doi: 10.1111/pan.12266. Epub 2013 Sep 25.
61. Pediatrics - Carr EA, Wilmoth ML, Eliades AB, Baker PJ, Shelestak D, Heisroth KL, Stoner KH (Akron Children's Hospital). Comparison of Temporal Artery to Rectal Temperature Measurements in Children Up to 24 Months, *Journal of Pediatric Nursing*, In Press, [Epub ahead of print], Jan 25, 2010.

62. Pediatrics - Health Canada (2017). Summary Safety Review - Ear and Forehead (contact) Infrared Thermometers (various brands) - Assessing the potential risk of inaccuracy in children under 2 years old. <https://www.canada.ca/en/health-canada/services/drugs-health-products/medefect-canada/safety-reviews/summary-safety-review-forehead-contact-infrared-thermometers-various-brands-assessing-potential-risk-inaccuracy-children-under-2-years-old.html>
63. Pediatrics - Yang WC, Kuo HT, et al. Tympanic temperature versus temporal temperature in patients with pyrexia and chills. *Medicine (Baltimore)*. 2016 Nov;95(44):e5267.
64. Pediatrics - Paul IM, Sturgis SA, Yang C, Engle L, Watts H, Berlin CM Jr (Penn State College of Medicine). Efficacy of standard doses of Ibuprofen alone, alternating, and combined with acetaminophen for the treatment of febrile children. *Clin Ther*. 2010 Dec;32(14):2433-40.
65. Pediatrics - Callanan D (Christus Santa Rosa Children's Hospital). Detecting fever in young infants: reliability of perceived, pacifier, and temporal artery temperatures in infants younger than 3 months of age. *Pediatr Emerg Care*. 2003 Aug;19(4):240-3.
66. Pediatrics - Allegaert K, Casteels K, van Gorp I, Bogaert G. Tympanic, infrared skin, and temporal artery scan thermometers compared with rectal measurement in children: a real-life assessment. *Curr Ther Res Clin Exp*. 2014 May 8;76:34-8. doi: 10.1016/j.curtheres.2013.11.005. eCollection 2014.
67. Pediatrics - Al-Mukhaizeem F, Allen U, Komar L, et al (University of Toronto/Hospital for Sick Children). Validation of the temporal artery thermometry by its comparison with the esophageal method in children. *Pediatric Academic Societies Annual Meeting*, May 3-6, 2003, Seattle, WA
68. Pediatrics - Al-Mukhaizeem F, Allen U, Komar L, et al (University of Toronto/Hospital for Sick Children). Comparison of temporal artery, rectal and esophageal core temperatures in children: Results of a pilot study. *Journal of Pediatric and Child Health*, Vol 9, No 7, pp 461-465, 2004
69. Pediatrics - Bahorski J, Repasky T, Ranner D, Fields A, Jackson M, Moultry L, Pierce K, Sandell M (Tallahassee Memorial Healthcare). Temperature measurement in pediatrics: a comparison of the rectal method versus the temporal artery method. In Press, Corrected Proof, Available online 24 February 2011, *Journal of Pediatric Nursing* (2011).
70. Pediatrics - Siberry GK, Diener-West M, Schappell E, Karron RA (Department of Pediatrics, School of Medicine, The Johns Hopkins University). Comparison of temple temperatures with rectal temperatures in children under two years of age. *Clinical Pediatrics*, pp 405-414, July/August 2002.
71. Pediatrics - Roy S, Powell K, Gerson LW (Akron Children's Hospital). Temporal artery temperature measurements in healthy infants, children, and adolescents. *Clinical Pediatrics*, pp 433-437, June 2003.
72. Peri-operative - Harper CM (Royal Sussex County Hospital Brighton, UK). The need for an accurate noninvasive thermometer. *Anesth Analg*. 2009 Jul;109(1):288; author reply 288-9.
73. Pharmacy - Canales AE (Texas Tech University Health Sciences Center). OTC device: temporal scanner TAT-2000C. *J Am Pharm Assoc (Wash DC)*. 2007 Jan-Feb;47(1):112.
74. School Nurse - Pappas M. Understanding the different methods for taking a temperature. *NASN School Nurse* 2012 27: 254 originally published online 5 July 2012.

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