

TAT-5000S Externally Powered Models

The following TAT-5000S externally powered models are authorized to be used with the monitors in the chart below.

Model	Exergen p/n	Authorized Monitors
TAT-5000S-RS232-PH	124265	Philips EarlyVue VS30 Philips IntelliVue MP 5SC
TAT-5000S-232EP	124229-AF-EP	Philips SureSigns VS4

The TAT-5000S models listed above are powered by an external 5V source provided by the patient/vital signs monitors listed above. Only the monitors listed are qualified to be used by the TAT-5000S models listed above.